- EN:** Instructions for Use
Ankura™ AAA and Cuff Stent Graft System
- BG:** Инструкции за употреба
стент-графт системи Ankura™: AAA и с маншон
- CS:** Návod k použití
systému Ankura™ AAA a Cuff stent-graft
- DA:** Brugsanvisning
Ankura™ AAA og Cuff-stenttransplantatsystemer
- DE:** Gebrauchsanweisung
Ankura™ AAA und Cuff-Stentgraftsyste
- EL:** Οδηγίες χρήσης
Συστήματα ενδομοσχευμάτων Ankura™ AAA και Cuff
- ES:** Instrucciones de uso
Sistemas de endoprótesis AAA y Cuff Ankura™
- FR:** Notice d'utilisation
Systèmes d'endoprothèse Ankura™ AAA et Cuff
- HR:** Upute za uporabu
Ankura™ AAA i Cuff Stent Graft sustavi
- IT:** Istruzioni per l'uso
Ankura™ Sistemi di endoprotesi AAA e Cuff
- NL:** Gebruiksaanwijzing
Ankura™ AAA en Cuffstent transplantaatsystemen
- PL:** Instrukcja obsługi
Systemy stentgraftu AAA i Cuff Ankura™
- PT:** Instruções de utilização
Ankura™ Sistemas de Endoprótese AAA e Cuff
- RO:** Instrucțiuni de utilizare
Sisteme de greță stent AAA și Cuff Ankura™
- SV:** Bruksanvisning
Stentgraftsystemet Ankura™ AAA och Cuff
- TR :** Kullanım Talimatları
Ankura™ AAA ve Cuff Stent Greft Sistemleri

	<p>Date of manufacture Дата на производство Datum výroby Produktionsdato Herstellungsdatum Ημερομηνία κατασκευής Fecha de fabricación Date de fabrication Datum proizvodnje Data di fabbricazione Productiedatum Data produkcji Data de fabrico Data fabricației Tillverkningsdatum Üretim tarihi</p>
	<p>Use-by Date Срок на годност Datum spotřeby Sidste anvendelsesdato Verfallsdatum Ημερομηνία λήξης Fecha de caducidad Date de péremption Rok upotrebe Data di scadenza Uiterste gebruiksdatum Termin ważności Prazo de validade A se utiliza înainte de Används före-datum Son Kullanma Tarihi</p>
	<p>Serial number Серийн номер Sériové číslo Serienummer Seriennummer Αριθμός σειράς Número de serie Numéro de série Serijski broj Numero di serie Serienummer Numer seryjny Número de série Număr de serie Serienummer Seri numarası</p>
	<p>Double sterile barrier system Система с двойна стерилна бариера Systém dvojité sterilní bariéry Dobbelt sterilt barrieresystem System mit doppelter Sterilbarriere Σύστημα διπλού αποστειρωμένου φραγμού Sistema con doble barrera estéril Système à double barrière stérile Dvostruki sterilni sustav barijere Sistema con doppia barriera sterile Dubbel steriel barrieresysteem System podwójnej bariery sterylnej Sistema de barreira estéril dupla Sistem de barieră sterilă dublă Dubbla sterila barriärsystem Çift steril bariyer sistemi</p>
	<p>Sterilized using ethylene oxide Стерилизовано с этиленов оксид Sterilizováno za použití ethylenoxidu Steriliseret med ethylenoxid Mit Ethylenoxid sterilisiert Αποστειρωμένο με αιθυλενοξειδίο Esterilizado con óxido de etileno Stérilisé à l'oxyde d'éthylène Sterilizirano etilen-oksidom Sterilizzato con ossido di etilene Gesteriliseerd met ethyleenoxide Wyjałowiono tlenkiem etylenu Esterilizado utilizando óxido de etileno Sterilizat cu oxid de etilenă Steriliserad med etylen Etilen oksitle sterilize edilmiştir</p>
	<p>batch code Код на партида Kód šarže Batchkode Chargennummer Κωδικός παρτίδας Código de lote Code de lot Šifra lota Codice del lotto Partijcode Kod partii Código do lote Cod lot satskod parti kodu</p>
	<p>Do not reuse За еднократна употреба Neroužívejte opakovaně Må ikke genbruges Nicht wiederverwenden Να μην επαναχρησιμοποιείται No reutilizar Ne pas réutiliser Nije za višekratnu uporabu Non riutilizzare Niet opnieuw gebruiken Nie używać ponownie Não reutilizar A nu se reutiliza Återanvänd inte Tekrar kullanmayın</p>

Consult Instructions for Use | Видте инструкциите за употреба | Přečtěte si návod k použití | Se brugsanvisningen | Gebrauchsanweisung beachten | Συμβουλευθείτε τις οδηγίες χρήσης | Consultar las instrucciones de uso | Lire la notice d'utilisation | Proučite Upute za uporabu | Consultare le istruzioni per l'uso | Raadpleeg de gebruiksaanwijzing | Patrz Instrukcja stosowania | Consultar as instruções de utilização | A se consulta instrucțiunile de utilizare | Se bruksanvisning | Kullanma Kılavuzuna Başvurun

Keep away from the sunlight | Да се държи далече от слънчева светлина | Uchovávejte mimo dosah slunečního záření | Må ikke udsættes for sollys | Von Sonnenlicht fernhalten | Μακριά από το ηλιακό φως | Mantener alejado de la luz solar | Conserver à l'abri de la lumière | Čuvajte dalje od sunčeve svjetlosti | Tenere al riparo dalla luce solare | Uit de buurt van zonlicht bewaren | Chronić przed promieniowaniem słonecznym | Manter afastado da luz solar | A se păstra ferit de lumina directă a soarelui | Utsätt inte för solljus | Güneş ışığından uzak tutun

Keep dry | Да се държи сухо | Udržujte v suchu | Skal holdes tør | Trocken halten | Να διατηρείται στεγνό | Mantener seco | Conserver au sec | Čuvajte na suhom mjestu | Conservare in luogo asciutto | Droog bewaren | Przechowywać w suchym miejscu | Manter seco | A se păstra ferit de umezeală | Förvaras torrt | Kuru ortamda saklayın

Catalogue number | Каталоген номер | Katalogové číslo | Katalognummer | Katalognummer | Αριθμός καταλόγου | Número de referencia | Référence | Kataloški broj | Numero di catalogo | Catalogusnummer | Numer katalogowy | Número do catálogo | Număr catalog | Artikelnummer | Katalog numarası

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
Read all instructions carefully. Failure to follow the instructions, warnings and precautions may lead to serious complications or injury to the patients.

1. Device Description


● Ankura™ Stent Graft System

- ◇ The Ankura™ AAA Stent Graft System includes:
 - ◆ The Ankura™ AAA Stent Grafts
 - ◆ The Ankura™ Delivery System
- ◇ The Ankura™ Cuff Stent Graft System includes:
 - ◆ The Ankura™ Cuff Stent Grafts
 - ◆ The Ankura™ Delivery System
- ◇ The Ankura™ Stent Graft is pre-loaded into the Delivery System. The loaded delivery system is inserted endoluminally via the femoral or iliac artery and tracked through the patient's vasculature to deliver the stent graft to the target site.

● Ankura™ Stent Grafts

(a) Ankura™ AAA Stent Graft

(b) Ankura™ Cuff Stent Graft

Fig. 1 Ankura™ AAA Stent Graft and Cuff Stent Graft


Table 1 Stent Graft Materials

Component	Material	Concentration in wt. %
Support Spring	Nitinol	49-67
Stainless steel sleeve	316LVM	7-17
Graft	PTFE	20-40
Marker	Platinum-Iridium	0.5-6


Suture	Polypropylene	<1
--------	---------------	----

✧ The stent grafts used in the Ankura™ AAA Stent Graft System are modular and comprises of three stent graft components, one Ankura™ AAA Stent Graft and two Ankura™ Cuff Stent Grafts. The modular system consists of one AAA Stent Graft and two Cuff Stent Grafts (Fig.1). Each component is introduced separately into the patient's vascular system. The Ankura™ Stent Graft is a tube composed of ePTFE film supported by a metal mesh called a stent. The markers are fixed on the stent. The graft which is ultrathin artificial vascular material fixed to the stent without suture, can isolate aneurysm from blood, and it is stronger than the weakened artery and it allows the blood to pass through it without pushing on the bulge. The stent can provide stable and expansionary force necessary to open the graft.

● Ankura™ Delivery System

(a) Ankura™ Delivery System (AAA Stent Graft)

(b) Ankura™ Delivery System (Cuff Stent Graft)


1	Tapered Tip	2	Markband	3	Proximal Anchor	4	Stent Graft
5	Inner Tube	6	Immobile Pole	7	Sheath	8	Front Grip
9	Slider Grip	10	Trigger	11	Screw Gear	12	Haemostatic Vale
13	Rear Grip	14	Proximal Releaser	15	Secure Buckle	16	Luer Connector
17	Front Cover	18	Rear Cover				

Fig.2 Ankura™ Delivery Systems


✧ The Ankura™ Stent Graft is pre-loaded inside the Ankura™ Delivery System. The Ankura™ Delivery System facilitates the placement of the stent graft via the arterial vasculature (e.g., femoral arteries). Under the fluoroscopic guidance, the Ankura™ Delivery System is properly positioned within the patient's vasculature and the stent graft is deployed from the delivery system. The AAA Stent Graft with Delivery System is shown in Figure 2(a), and the Cuff Stent Graft with Delivery System is shown in Figure 2(b).

● Device specifications

✧ The specification of Ankura™ AAA and Cuff Stent Graft System are indicated as follow and the details are shown respectively in table 2 and Table 3.

(a) AAA Stent Graft System

(b) Cuff Stent Graft System

Table 2 Summary of Ankura™ AAA Stent Graft System

Proximal Diameter (D1/mm)	Distal Diameter (D2/mm)	Distal Diameter (D3/mm)	Covered Length (L1/mm)	Covered Length (L2/mm)	Covered Length (L3/mm)	Delivery System	Recommended Vessel Diameter (mm)
20	10	10	100,110,120,130,140,160	50	30	20F	16-18
22			100,110,120,130,140,160			20F	18-20
24	12	12	100,110,120,130,140,160			20F	20-22
26			100,110,120,130,140,160			20F	21-24
28			100,110,120,130,140,160			20F	23-25
30			100,110,120,130,140,160			22F	25-27
32			100,110,120,130,140,160			22F	26-29
34			100,110,120,130,140,160			22F	28-31
36			100,110,120,130,140,160			22F	30-33

Note: the attached table is only the recommended reference size, and the final selection of stent graft specifications is subject to the doctor's judgment based on the actual situation of the patient.

Table 3 Summary of Ankura™ Cuff Stent Graft System

AAA Specification	Proximal Diameter (D1/mm)	Distal Diameter (D2/mm)	Covered Length (L1/mm)	Delivery System	Distal Recommended Vessel Diameter (mm)
20~22	12	10	40,60,80,90,100,120	18F	8-9
		12	40,60,80,90,100,120	18F	10-11
		14	40,60,80,90,100,120	18F	11-12
		16	40,60,80,90,100,120	18F	13-14
		18	40,60,80,90,100,120	18F	15-16
24~36	14	10	40,60,80,90,100,120	18F	8-9
		12	40,60,80,90,100,120	18F	10-11
		14	40,60,80,90,100,120	18F	11-12
		16	40,60,80,90,100,120	18F	13-14
		18	40,60,80,90,100,120	18F	15-16
		20	40,60,80,90,100,120	18F	16-18
		22	40,60,80,90,100,120	18F	18-20

Note: the attached table is only the recommended reference size, and

the final selection of stent graft specifications is subject to the doctor's judgment based on the actual situation of the patient.

Physicians should strictly adhere to indications and contraindications.

2. Indications for Use

- Ankura™ AAA Stent Graft System is indicated for endovascular repair of patients with abdominal aortic aneurysms, requiring:
 - ◇ Adequate iliac/femoral access vessel that is compatible with the required delivery system.
 - ◇ Non-aneurysmal proximal aortic neck length ≥ 15 mm.
 - ◇ Non-aneurysmal proximal aortic neck diameter of 18–32mm.
 - ◇ Proximal aortic neck angulation $\leq 60^\circ$.
 - ◇ Distal iliac artery anchorage zone ≥ 15 mm.
 - ◇ Distal iliac artery diameter of 8–20mm.
 - ◇ Morphology suitable for endovascular repair.
- Additional considerations for patient selection include:
 - ◇ Patient's age and life expectancy;
 - ◇ Co-morbidities (e.g., cardiac, pulmonary or renal insufficiency prior to surgery, morbid obesity);
 - ◇ Toleration to the general, regional or local anesthesia;
 - ◇ The final decision is at the discretion of the physician and patient.
- Device Selection
 - ◇ The size of the Ankura™ Stent Graft should be appropriate to fit the patient's aortic anatomy. Lifetech recommends an oversized Ankura™ Stent Graft by 10%-20%.

3. Contraindications

- The Ankura™ AAA Stent Graft is contraindicated in:
 - ◇ Patients with acute systemic infection;
 - ◇ Patients with mesenteric blood flow mainly supplied by the inferior mesenteric artery;
 - ◇ Patients with accessory renal artery original from abdominal artery;
 - ◇ Patients who have allergic reaction to the device;
 - ◇ Patients who are not suitable for endovascular repair in vascular morphology;
 - ◇ Patients who can not tolerate contrast agents due to renal insufficiency;
 - ◇ Patients who are allergic to contrast agents;
 - ◇ Aneurysms neck with thrombus.
 - ◇ Non-aneurysmal proximal aortic neck length < 15 mm.
 - ◇ Non-aneurysmal proximal aortic neck diameter < 18 mm or > 32 mm.
 - ◇ Proximal aortic neck angulation $> 60^\circ$.
 - ◇ Distal iliac artery anchorage zone < 15 mm.
 - ◇ Distal iliac artery diameter < 8 mm or > 20 mm.

4. Pre-implantation Preparation

Read the instructions carefully before using the device. The following instructions embody a basic guideline for device placement. Variations in the following procedures may be necessary. These instructions are intended to help guide the physician rather than replacing the physician judgment.

- Materials required
 - ◇ 0.035"/260cm Super Stiff Guidewire. For example, Lifetech's Super Stiff Guidewire, Amplatz Super Stiff Guidewire or comparable.
 - ◇ Multi-plane fluoroscopic imaging system;
 - ◇ Angiographic sheath and marker sheath;
 - ◇ Contrast agent;
 - ◇ High pressure syringe;
 - ◇ Heparinized saline solution;
 - ◇ Any other devices required by the physician for the procedure (including devices used for femoral artery dissection).
- A multidisciplinary team that has combined procedural experience with:
 - ◇ Local anesthesia or general anesthesia;
 - ◇ Femoral artery exposure, arteriotomy and repair;
 - ◇ Percutaneous access and closure techniques;
 - ◇ Fluoroscopic and angiographic image interpretation;
 - ◇ Endovascular stent placement;
 - ◇ Appropriate use of radiographic contrast agent;
 - ◇ Techniques to minimize radiation exposure;

- ✧ Expertise in necessary patient follow-up modalities;
- ✧ Experience in traditional vascular surgery for any possibility of failure in EVAR.
- **Products preparation**
 - ✧ The Stent Graft diameter should be selected by the vessel inner diameter, and not outer diameter. The improper size of the device may result in endoleak or stent-graft migration. Generally, the diameter of Stent Graft selected should be 1.1 to 1.2 times larger than the landing zone vessel inner diameter.
 - ✧ Inspect the device and the package to verify that no damage has occurred as a result of shipping. Do not use the device if any damage of the device or the sterilization barrier has been observed. If damage has occurred, return it to Lifetech. Prior to use the device, re-review the specification conforming to what the patient requires.
- **Patient preparation**
 - ✧ Measure proximal aortic neck length and vessel inner diameter.
 - ✧ Measure the angle of the proximal aortic neck relative to the long axis of aneurysm.
 - ✧ Evaluate the quality of aortic neck.
 - ✧ Assess the morphology of the vessel compatible with endovascular aortic repair.
 - ✧ Undergo any necessary common examination, and monitoring of vital signs.
 - ✧ Refer to institutional protocols relating to anesthesia, anticoagulation and monitoring of vital signs.
 - ✧ Position patient on imaging table allowing fluoroscopic visualization from the aortic arch to the femoral artery bifurcations.
 - ✧ Adequate iliac/femoral arteries access.
 - ✧ Expose both common femoral arteries access using standard surgical technique.
 - ✧ Systemic anticoagulation should be used during the deployment procedure based on hospital and physician preferred protocol. If heparin is contraindicated, an alternative anticoagulant should be considered.

5. Implantation Procedure

- **Preparation of AAA and Cuff delivery system before implantation**
 - ✧ Attach the syringe with heparinized saline to the luer hub and hemostasis valve and flush separately to evacuate air from the inner tube and the outer sheath. (Fig.3)


Fig 3 Flush the inner tube and the outer sheath

- ✧ Monitor important vital signs.
- ✧ Administer heparin to adjust ACT(activated clotting time) to be 200 seconds. Inspect flush solution. Sheath and guidewire should be flushed before every insertion.
- ✧ Use gauze pads with saline solution to wipe the outer sheath and activate the hydrophilic coating. (Fig. 4)


Fig. 4 Wipe the sheath with sterile heparinized normal saline

- ◇ Expose bilateral femoral arteries by standard surgery.
- ◇ Insert one 0.035"/260cm super stiff guidewire in the femoral artery, and advance it until above the renal arteries.
- ◇ Insert marker sheath, and double check the stent specification by radiography.
- Deploy the stent graft
 - ◇ Advance the delivery system to the abdominal aneurysm along the guidewire, and always maintain the hemostasis valve towards the other side of the patient. (Fig. 5)


Fig. 5 Introduction of the AAA Stent Graft with Delivery System

Note: The Short branch of the AAA Stent Graft must be in the same side with the hemostasis valve of the delivery system when the AAA Stent Graft is preloaded into the delivery system.

- ◇ Inject contrast agents into the abdominal aorta; mark the position of the target location, either on the imaging screen or with radiopaque marker or angiographic ruler on the patient's body. Adjust the position of the Stent Graft through the proximal radiopaque markers so that the top edge of the graft is just below the lowest renal artery.
- ◇ Ensure that the proximal radiopaque "O-shaped" marker is on the same side with the "V-shaped" marker on the long branch of the AAA stent graft: Ensure that the proximal radiopaque "8-shaped" marker is on the same side with the "e-shaped" marker, the "O-shaped" marker and the "✓-shaped" marker on the short branch of the AAA stent graft.
- ◇ Ensure that the distal portion of the short branch is above the aortic bifurcation and within the aneurysmal sac.
- ◇ Slowly deploy the stent graft. Hold the front grip of the delivery system stationary, then rotate the slider grip counterclockwise to slowly deploy the stent graft until the short branch is deployed. (Fig. 6)


Fig. 6 Deploy the short branch of AAA Stent Graft

- Proximal end release

- ◇ Keep the delivery system stationary. Remove the Secure Buckle from the grip of the delivery system. (Fig. 7)


Fig. 7 Remove the Secure Buckle

- ◇ Pull the Proximal Releaser to release the proximal spring of the stent graft. (Fig. 8)



Fig. 8 Pull the Proximal Releaser

- Insert another guidewire through the other side of the femoral artery into the short branch of placed AAA Stent Graft. (Fig. 9)


Fig. 9 Introduction of another guidewire

- ◇ Check whether there is a gap between tip head and sheath of Cuff delivery system, if necessary, loosen the rear cover, gently pull back sheath core until there is no gap, and then lock the rear cover.
- ◇ Check whether the rear cover of Cuff delivery system is loose, it is necessary to tighten before subsequent operation.
- After confirming the rear cover of Cuff delivery system is locked, insert the Cuff stent graft with delivery system over guidewire into the short branch of the AAA stent graft, and always maintain the haemostatic valve towards the inside of the patient. As show in (Fig. 10). Align the radiopaque marker on the Cuff stent with the radiopaque marker on the short branch of the AAA Stent Graft, with at least two waves overlapped between them. Ensure that the radiopaque "O-shaped" marker on the Cuff stent graft is almost as high as the "e-shaped" marker on the short branch of the AAA Stent Graft.


Fig. 10 Insert the Cuff Stent Graft with Delivery System

- Hold the front grip stationary, unlock the front cover anticlockwise (Fig. 11).


Fig. 11 Unlock the front cover anticlockwise

- Hold the rear grip stationary with one hand and pull back the front grip with the other to deploy the extension branch Cuff under monitoring (Fig. 12). After the stent graft is deployed, hold the front grip stationary, and pull back the rear grip to make the tip back to the distal end of the out sheath. Then, withdraw the delivery system and guidewire. (Fig. 13)


Fig. 12 Delivery System Status of Pulling Back Front Grip


Fig. 13 Deploy the Cuff Stent Graft in the short branch of AAA Stent Graft

- Rapidly deploy the long branch of the AAA Stent Graft. Hold the front grip of the delivery system stationary, then unlock the trigger and pull the slider grip along the screw gear to fully deploy the stent graft. (Fig. 14)
- Hold the slider grip of the delivery system stationary, unlock the trigger and pull back the front grip to make the tip back to the distal end of the out sheath. Then remove the delivery system.


Fig. 14 Deploy the long branch of the AAA Stent Graft

- Insert another Cuff stent graft system over guidewire into the long branch of the AAA Stent Graft, and always maintain the haemostatic valve towards the inside of the patient. As show in Fig. 15. Align the radiopaque marker on the Cuff Stent Graft with the radiopaque marker on the long branch of the AAA Stent Graft, with at least two waves overlapped between them. Ensure that the radiopaque "O-shaped" marker on the Cuff stent graft is almost as high as the "o-shaped" marker on the long branch of the AAA Stent Graft. (Fig. 15)


Fig. 15 Insert another Cuff Stent Graft

- Hold the front grip stationary, unlock the front cover anticlockwise. (Fig. 11) Hold the rear grip stationary with one hand and pull back the front grip with the other to deploy the extension branch Cuff under monitoring (Fig. 12). After the stent is deployed, hold the front grip stationary, and pull back the rear grip to make the tip back to the distal end of the out sheath. Then, withdraw the delivery system and guidewire. (Fig. 16)


Fig. 16 Deploy the Cuff Stent Graft in the long branch of AAA Stent Graft

- Handle Disassembly Technique for Partial Stent Graft Deployment
In the unlikely event of delivery system failure and concomitant partial stent graft deployment, a "handle disassembly" technique may permit the successful deployment of the stent graft.
Step 1 - Pull back the trigger and fully retract the slider. Note: Since the graft cover is severed, the slider can be retracted without further deploying the stent graft.
Step 2 - Stabilize the delivery system.
Step 3 - Insert the tips of a pair of hemostats into gap at the end handle on the slider grip. Then, slide the hemostatic forward to split the slider grip into left and right halves.
Step 4 - Remove the left and right halves of the slider grip from the screw gear to expose the sheath joint.
Step 5 - Hold the front handle with one hand and hold the sheath joint with the other hand and pull back until the stent graft is fully deployed.
Step 6 - Remove the delivery system by gripping the screw gear

and withdrawing from the patient.

6. Potential Adverse Events

Adverse events associated with use of the Ankura™ AAA Stent Graft System include, but are not limited to the following:

- Allergic reaction
- Aneurysms enlargement
- Aneurysms/vessel rupture
- Anesthetic complications
- Arrhythmia
- Arterial or venous thrombosis (including thrombosis and embolization)
- Arterial damage
- Arterial perforation
- Arteriovenous fistula or pseudo aneurysm
- Bleeding/Hematoma
- Bowel complications
- Cardiac failure/Infarction
- Coagulopathy
- Conversion to open surgery
- Death
- Edema
- Erectile dysfunction
- Embolism
- Hemorrhage/bleeding
- Infection
- Intercostal pain
- Neurological damage, local or systemic (e.g., stroke, paraplegia, paraparesis)
- Paraparesis
- Paralysis
- Paresthesia
- Pulmonary complications
- Renal function insufficiency/kidney failure
- Shock
- Tissue necrosis
- Wound complications and subsequent attendant problems
- Device complications due to:
 - ◇ Insertion and removal difficulties
 - ◇ Release difficulties
 - ◇ Failure to deploy
 - ◇ Inaccurate placement
 - ◇ Endoleak
 - ◇ Stent migration
 - ◇ Stent fracture
 - ◇ Rupture of graft material
 - ◇ Erosion
 - ◇ Stent thrombosis or embolization
 - ◇ Stent graft expansion
 - ◇ Stent graft twisting or kinking

Note:

For any serious incident, the user or patient should report it to Lifetech at quality@lifetechmed.com or the local sale representative of Lifetech and your local Competent Authority (CA).

7. Summary of Clinical Data

There were two complete clinical studies of Ankura™ AAA/ AUI stent graft system, one in China and the other in Turkey. Both studies aimed to evaluate the clinical safety and efficacy of the Ankura™ Stent Graft Systems, and they met the predefined primary endpoints of treatment success with few adverse events.

Study in China

107 patients implanted Ankura™ stent graft systems (TAA-62 cases, AAA-30 cases, AUI-15 cases) were enrolled in the study in China. The stent grafts were successfully implanted in all 107 patients (100%). 91 of 106 patients were followed up after the procedure, and the mean follow-up time was 218 (range 30 – 600) days. Treatment success rate of aortic aneurysm at 6-month post-implant reached 93.4% (85/91). 85 of 91 patients evaluated for safety events during the follow up were free of a safety event (93.4%).

Follow-up data showed that:

One out of 107 patients (0.9%) died in-hospital (0.9%), but this death was unrelated to the device or procedure.

One out of 107 patients (0.9%) had artery stenosis.

One out of 107 patients (0.9%) experienced paraplegia.

One out of 107 patients (0.9%) had renal failure.

Three out of 107 patients (2.8%) experienced early endoleaks.

No patient experienced stent displacement, stent fracture, or aortic rupture.

Study in Turkey

The study in Turkey recruited 100 patients and the immediate treatment success rate was 100%. These patients were evaluated for safety events through one year.

One-year patient follow-up showed that:

Five patients (5%) experienced type 2 endoleak.

Two patients (2%) had type 1b endoleak due to native iliac artery dilatation.

Two patients (2%) experienced secondary interventions.

Two patients (2%) died due to cardiac and neurological causes.

No patient experienced device migration, device deficiency, device replacement, stent kink or fracture, limb occlusion or conversion to open surgery.

8. Storage and Transport Requirements

- The Ankura™ AAA Stent Graft and Cuff Stent Graft are pre-loaded in delivery system, and stored in sterile package.
- This device should be stored in a clean, well ventilated room without corrosive gas and stored at room temperature.
- During transportation, it is strictly prohibited to squeeze, get wet and expose to the sun, and handle with care.
- No exposure to organic solvent.

9. Warnings

● General

- ◇ Read all instructions carefully. Failure to properly follow the instructions, warnings and precautions may lead to serious consequences or injury to the patient.
- ◇ The Ankura™ AAA Stent Graft and Cuff Stent Graft cannot be used after resterilization.
- ◇ If any sign of damage is observed prior to use, do not use the product and return it to Lifetech.
- ◇ This device should only be used by physicians and teams trained in vascular interventional techniques and the use of this device.
- ◇ The long term performance of endovascular grafts has not yet been confirmed. All patients should be advised that endovascular treatment requires lifelong, regular follow-up to assess their health and the performance of their endovascular graft (e.g., endoleak, migration, distortion, aneurysms enlarging etc).
- ◇ Clinical studies show that improper patient selection may result in poor treatment effect. In addition, unsuccessful endovascular aortic repair which requires open surgery may result in an increased risk of complications related to both procedures.
- ◇ Patient's condition needs to be evaluated before implantation on the following aspects, cardiac failure/infarction, pulmonary complications, edema, and renal function insufficiency.
- ◇ After the device placement, patients should be regularly monitored for the state of the device, including endoleak, migration, change of the Stent Graft structure or position, aneurysms enlarging etc. At a minimum, annual X-ray, contrast or non-contrast CT imaging is required.
- ◇ Patients with specific clinical findings, e.g, enlarging aneurysms or unexpected diseases in the position of the endovascular graft should undergo secondary interventions or open surgical treatment. An increase in aneurysms size and/or persistent endoleak may lead to aneurysms rupture.
- ◇ Patients with infection after the device placement should be treated with oral antibiotics.
- ◇ A vascular surgery team should always be available during the implantation, or re-intervention procedure in which situation that open surgical repair is needed.
- ◇ Do not re-use or re-sterilize the Ankura™ AAA Stent Graft System and the Ankura™ Cuff Stent Graft System. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or cause device failure which, in

turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also cause a risk of contamination of the device and/or patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

- Patients selection, treatment and follow-up
 - ◇ The safety and effectiveness of the Ankura™AAA Stent Graft has not yet been evaluated in following populations:
 - ◆ Traumatic aortic injury
 - ◆ Uncorrectable coagulopathy
 - ◆ Hereditary connective tissue disease (e.g., Marfan's or Ehlers-Danlos Syndromes)
 - ◆ Concomitant thoracic aortic or thoraco abdominal aneurysms
 - ◆ Patients with active systemic infections
 - ◆ Pregnant or nursing females
 - ◆ Morbid obesity
 - ◆ Less than 18 years of age
 - ◆ Patient's life expectancy less than one year
 - ◇ Access vessel diameter (measured inner wall to inner wall) and morphology (minimal tortuosity, occlusive disease and/or calcification) should be compatible with vascular access techniques and the profile of the delivery system. Too tortuous vasculature may increase the risk of failure in stent graft deployment or proximal end release.
 - ◇ The Ankura™AAA Stent Graft System is not recommended in patients who weigh greater than 350 Pounds (150kg) or cannot undergo accurate fluoroscopy examination due to obesity.
 - ◇ Patients with a systemic infection may be at increased risk of endovascular stent-graft infection.
 - ◇ Patients with uncorrectable coagulopathy may have an increased risk of Type II endoleak or bleeding complications.
 - ◇ Patients need to measure blood pressure every day after surgery and make a clinical record. High blood pressure may lead to stent expansion, vascular complications, and other complications. Blood pressure should be monitored in time for the medical staff to grasp the postoperative status of the patient. At the same time, the doctor needs to remind the patient to control blood pressure in the normal range.
 - ◇ The product should not be used in patients who are allergic to the materials listed out in Table 1.
- Implant Procedure
 - ◇ Systemic anticoagulation should be used during the procedure based on hospital and physician preferred protocol. If heparin is contraindicated, an alternative anticoagulant should be considered.
 - ◇ Minimize handling of the constrained endoprosthesis during preparation and insertion to decrease the risk of endoprosthesis contamination and infection.
 - ◇ Maintain the guide wire position during delivery system insertion.
 - ◇ Do not bend or kink the delivery system. Doing so may cause damage to the delivery system and the stent graft.
 - ◇ Always use fluoroscopy for guidance, delivery and observation of the device within the vasculature.
 - ◇ The use of the Ankura™ Stent Graft System requires administration of intravascular contrast. Patients with preexisting renal insufficiency may have an increased risk of renal failure postoperatively. Care should be taken to limit the amount of contrast agent used during the procedure.
 - ◇ To avoid any twist in the endovascular stent-graft, during any rotation of the delivery system, rotate all of the components of the system together carefully.
 - ◇ Inaccurate placement of the stent graft within the vessel may result in increased risk of endoleak, migration or embolization.
 - ◇ Inadequate fixation of the stent graft may result in increased risk of migration of the stent graft, and incorrect deployment and migration of the stent graft requires surgical intervention.
 - ◇ Do not continue advancing if resistance is felt during advancement of the guide wire or delivery system. Stop and assess the cause of resistance.
 - ◇ Unless medically indicated, do not deploy the stent graft in a

location that will occlude arteries necessary to supply blood flow to organs or extremities.

- **MRI Information**

Nonclinical testing has demonstrated that the Ankura™ Stent Graft is MR Conditional. A patient with this device can be scanned safely under the following conditions:

- ◇ Static magnetic field of 3.0 Tesla or less
- ◇ Spatial gradient field of ≤ 720 Gauss/cm
- ◇ Maximum whole-body-averaged specific absorption rate (SAR) of 2W/kg for 15 minutes of scanning

Based on non-clinical testing, the device was determined to produce a temperature rise of less than 2°C at a maximum whole-body-averaged specific absorption rate (SAR) of 2W/kg for 15 minutes of MR scanning. The maximum whole-body-averaged specific absorption rate (SAR) was derived by calculation and verified by calorimetry.

Note: MR image quality may be compromised if the area of interest is in the same area or relatively close to the position of the device. Therefore, it may be necessary to optimize MR imaging parameters for the presence of this implant. The image artifact extends approximately 20 mm outside the device lumen when scanned in non-clinical testing using the spin echo sequence and gradient echo sequence, respectively in a 3.0 T Scanner with a whole body coil.

10. Patient consultation information

The biggest benefits of treatment of the patient's diseased or injured aorta with Ankura™ stent graft are decreased chance of rupture and restoration of normal blood flow. If left untreated, aortic lesions can expand and rupture, resulting in bleeding inside the body, which is life-threatening.

The physician should consider, but not be limited to, the following when explaining this endovascular device and related procedures to the patient:

- The differences between endovascular repair and open surgical repair, include the risks of open surgical repair and endovascular repair.
- Pros and cons of endovascular repair and open surgery repair.
- Endovascular repair has the potential advantage of minimal trauma.
- Endovascular repair or open surgical repair of the lesion may be required in the future.
- Details of possible complications after implantation of the device.
- Regular follow-up after implantation of the device.

11. Follow Up

- Periods of follow up
 - ◇ Clinical follow up should be performed at 24 hours, 1, 3, 6, 12 and 24 months after the procedure. The follow up period can be adjusted at anytime according to the circumstances of each individual patient.
- Methods used during follow up
 - ◇ Haematological examination, X-ray, B-US, CTA, MRI, DSA etc.
- Follow up items:
 - ◆ Size of aneurysms;
 - ◆ Embolization;
 - ◆ Changing of aneurysms pulsatility;
 - ◆ Migration;
 - ◆ Endoleak;
 - ◆ Stent Graft distortion.
- ◇ In addition, patients with following specific findings during following up, endovascular intervention or surgical treatment should be considered:
 - ◆ Aneurysms enlarged 5mm and above;
 - ◆ Changing of aneurysms pulsatility with or without aneurysms enlarging or endoleak;
 - ◆ Continual endoleak with or without aneurysms enlarging;
 - ◆ Endoleak caused by Stent Graft migration.

12. Shelf Life

- The Ankura™ AAA Stent Graft System and Ankura™ Cuff Stent Graft System are sterilized with ethylene oxide. The shelf life is indicated on the label. Do not use the expired product.

13. Package and Label

- The Ankura™ AAA Stent Graft System and Ankura™ Cuff Stent Graft System are supplied sterile and packaged separately. Each Ankura™ Stent Graft is individually contained within a Delivery System.
- The Ankura™ AAA and Cuff Stent Graft with Delivery Systems are immobilized and protected in a PETG tray respectively, then sealed with two Tyvek1073B pouches, on which a primary label is stuck. The product is sterilized and put in a box with IFU, implant card, customer feedback form and other accompanying document. A label is stuck on the box.

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Lifetech Scientific (Shenzhen) Co., Ltd