

Self-expanding Peripheral Stent 635

EXIST

The Qmedics portfolio is based on a “back to the basic” concept, where each product has been carefully designed to be compliant with specific lesions and movability of the artery.


The Exist Flex provides a perfect balance between radial compression resistance and flexion deformation. Recommended to be used in complex and moderate calcified lesions mostly subjected to bending deformation.

The Exist Pull provides perfect balance between compression resistance and longitudinal conformability. Recommended to be used in highly calcified, challenging and difficult to cross lesions within arteries subjected to longitudinal deformations and radial compression

EXIST STENT 635

- Diameter 5 to 10 mm
- Stent length 30 to 200 mm
- Pull-back OTW delivery system
- Multifunctional braided catheter shaft
- Radiopaque tip
- Marker band on outer shaft for accurate stent deployment
- Usable length 85, 135 and 175 cm
- Strut thickness $\leq 140 \mu\text{m}$
- 2 tantal markers at each stent end
- 6 F introducer compatible
- 0.035" guide wire compatible
- Peak-to-Valley stent design, extremely flexible and kink resistant
- Intuitive delivery system
- EtO-sterilised

RADIAL RESISTANCE WITH AXIAL CONFORMABILITY



Benchmarked at Qmedics AG

AVAILABLE SIZES

■ 85/135/175

Ø/L	30	40	60	80	100	120	140	160	180	200
5.0	■	■	■	■	■	■	■	■	■	■
6.0	■	■	■	■	■	■	■	■	■	■
7.0	■	■	■	■	■	■	■	■	■	■
8.0	■	■	■	■	■	■	■	■	■	■
9.0	■	■	■	■	■	■				
10.0	■	■	■	■	■	■				



Meet us online!

Qmedics AG | Gewerbezentrums Arova | Winterthurerstrasse | CH-8247 Flurlingen
 Phone +41 (0)44 515 82 02 | E-Mail info@qmedics.ch | Website www.qmedics.ch

