


# SeQure™ Snare System



**Multiple Use, Precise Control**

# FEATURES & BENEFITS

Figure : 90°angle between Loop and Shaft



- Titanium Nitride Coated Tungsten Loop
- Flexible Nitinol shaft provides Kink Resistance and 1:1 Torque control
- 90°angle between Loop and Shaft makes manipulation and retrieval more simple
- Radiopaque Markband at distal tip of catheter allows excellent Fluoroscopic Visualization

## ORDERING INFORMATION

### • Snare System-General

Code	Loop Diameter	Shaft Length	Catheter Size	Catheter Length
	mm	cm	Fr	cm
LT- SG- 05	5	120	4	100
LT- SG- 10	10	120	4	100
LT- SG- 15	15	120	6	100
LT- SG- 20	20	120	6	100
LT- SG- 25	25	120	6	100
LT- SG- 30	30	120	6	100
LT- SG- 35	35	120	6	100

### • Snare System-Petite

Code	Loop Diameter	Shaft Length	Catheter Size	Catheter Length
	mm	cm	Fr	cm
LT- SP- 10	10	65	4	50
LT- SP- 25	25	65	6	50

### • Snare System-Micro

Code	Loop Diameter	Shaft Length	Catheter Size	Catheter Length
	mm	cm	Fr	cm
LT- SM- 02	2	200	3	175
LT- SM- 04	4	200	3	175
LT- SM- 07	7	200	3	175

For further inquiries, please contact your local sales representatives.

### Lifetech Scientific (Shenzhen) Co., Ltd.

Add: Floor 1-5, Cybio Electronic Building, Langshan 2nd Street, North Area of High-tech Park, Nanshan District, Shenzhen 518057, P. R. China

Tel: +86-755-86026250 Fax: +86-755-86026251

Email: marketing@lifetechmed.com

Web: www.lifetechmed.com

EU authorized representative:  
Life Tech Scientific (Europe) Coöperatief U.A.  
Oliemolenstraat 60, 6416CB Heerlen, the Netherlands